

Routine analysis of glucose, urea and lactate

We perform routine analysis of glucose urea and lactate for you. Your laboratory will be disburdened and you will gain free capacities for your research.

With our pool of well proven and high quality analytical devices we perform an accurate, sure and reliable determination of the analytes concentration of your samples.

Urea	
Instrument	Beckman BUN Analyzer 2
Min. concentration *	0,5 mmol/l
Max. concentration **	35,8 mmol/l
Accuracy	2,5 % at 17,9 mmol/l
Required sample volume ***	10 µl per sample
Glucose	
Instrument	Beckman Glucose Analyzer 2
Min. concentration *	0,6 mmol/l
Max. concentration **	24,9 mmol/l
Accuracy	2 % at 8,3 mmol/l
Required sample volume ***	10 µl per sample
Lactate	
Instrument	Biosen 5030 AutoCal
Min. concentration *	0,5 mmol/l
Max. concentration **	40,0 mmol/l
Accuracy	1,5 % at 12 mmol/l Lactat
Required sample volume ***	20 µl per sample

* Lower concentrations may be measured by means of an increased sample volume. This procedure requires individual tests.

** Higher concentrations may be measured by means of sample dilution. This procedure requires individual tests.

*** Volume for single determinations.

For your comfort we propose following procedure:

We offer you to test some of your samples in advance of your order to assure that your samples can be analyzed.

You will obtain the results as a measurement report (hardcopy and data) as a well founded basis for your further research and evaluation.



Contact and Ordering Information

Please feel so free to contact us if you have any question or if you need more information about our service. We like to help you with more information about routine analysis of glucose, urea and lactate.

To contact us or to request a quotation you may use the contact form on our homepage www.analysisio.com or write us an email to info@analysisio.com. Of course you may contact us also via phone or fax.

We look forward to plan with you the extend and procedure individually optimized to your demands..

You find an actual pricelist in the download area where you found this information.

Imprint

analysisio GmbH
Brandteichstrasse 20
D-17489 Greifswald, Germany
Tel.: +49 (0)3834 550263
Fax: +49 (0)3834 550222
E-mail: info@analysisio.com
General manager: Jens Pörksen
Commercial registry: Stralsund HRB 7139
Sales tax id. number: DE252491041

