

Laboratory consumables and supplies

BIOLOGIX lab consumables and supplies
Superior product quality - unbeatable prices



BIOLOGIX
RESEARCH COMPANY



BIOLOGIX laboratory supplies


BIOLOGIX Research Company is committed to the design and manufacturing of high quality scientific and biomedical laboratory supplies. In this catalogue you will find superior products at unbeatable prices:

- Centrifuge tubes page 2
- Microcentrifuge tubes page 3
- PCR tubes page 3
- Microtube racks page 4
- Cryogenic storage boxes page 4
- Cardbord freezer boxes page 5
- Pipet tips page 6
- Transfer pipets page 7
- Solution basins page 8
- Drosophila vials and bottles page 9
- Tri-cornered beakers page 10
- Petri dishes page 10
- Cell scrapers and lifters page 11
- Embedding cassettes page 11
- Slide storage boxes page 12
- Slide mailers page 12
- Contact and ordering information page 13

Contact

analysisio GmbH
Brandteichstrasse 20
D-17489 Greifswald

Fon: +49(0)3834 550263
Fax: +49(0)3834 550222
Mobile: +49(0)171 3406464
Email: info@analysisio.com



Centrifuge tubes



15 or 50 ml volume
 Sterilized with gamma-radiation, non-pyrogenic
 Polypropylene
 Plug-seal or flat top cap
 Conical bottom
 White, easy-to-read, printed graduations
 White writing area, for easy sample identification
 Maximum recommended RCF rating: 8400 x g
 0°C to 40°C
 Dimensions: 17 x 130 mm and 30 x 115 mm
 Packed as bulk or on a recyclable polyfoam rack



15 ml volume
 Polypropylene
 Autoclavable
 Molded graduation
 Conical bottom
 Non-sterile
 Press cap closures (caps must be ordered separately)
 Maximum recommended RCF rating: 3000 x g
 Urine analysis as well as general use
 Dimensions: 17 x 110mm
 Bulk pack of 500 tubes

Article number	Description
10-0151	15 ml, plug seal cap, sterile, 50 tubes/rack
10-0152	15 ml, plug seal cap, sterile, 50 tubes/bag
10-0501	50 ml, plug seal cap, sterile, 25 tubes/rack
10-0502	50 ml, plug seal cap, sterile, 50 tubes/bag
10-9151	15 ml, flat top cap, sterile, 50 tubes/rack
10-9152	15 ml, flat top cap, sterile, 50 tubes/bag
10-9501	50 ml, flat top cap, sterile, 25 tubes/rack
10-9502	50 ml, flat top cap, sterile, 50 tubes/bag
10-0122	15 ml, nonsterile, autoclavable, bulk pack of 500/bag
10-0123	Press cap closures for 10-0122, PE, bulk pack of 1000/bag

Contact

analysisio GmbH
 Brandteichstrasse 20
 D-17489 Greifswald

Fon: +49(0)3834 550263
 Fax: +49(0)3834 550222
 Mobile: +49(0)171 3406464
 Email: info@analysisio.com




Microcentrifuge tubes



RNase & DNase free
 Medical grade polypropylene
 Autoclavable
 Attached flat-top Lid
 RCF rating >13,000 x G
 Working temperature: 0 to 100C
 Freeze proof, boil Proof
 Molded graduations
 Shaded writing area
 Red, yellow, blue color or clear
 Bulk pack 500 tubes/bag,

PCR tubes



Pure polypropylene, RNase & DNase free
 0,2ml or 0,5 ml volume
 Thinwall tubes, conical Base
 Single tubes with attached caps or
 Attached in a row of 8 or 12 tubes with
 Dome-shaped caps, attached in a row of 8 or 12
 Packed in resealable bags
 RCF rating: 6000 x G
 Accelerate temperature transfer
 Improved sample stability
 Suitable for all common thermocyclers

Article number	Description
80-1500	Microcentrifuge tubes, 1,5 ml, clear bulk pack of 500/bag
80-1501	Microcentrifuge tubes, 1,5 ml, red color, bulk pack of 500/bag
80-1502	Microcentrifuge tubes, 1,5 ml, yellow color, bulk pack of 500/bag
80-1503	Microcentrifuge tubes, 1,5 ml, blue color, bulk pack of 500/bag
60-0005	PCR tube, attached cap, 0,5 ml, bulk of 1000 tubes/bag
60-0002	PCR tube, attached cap, 0,2 ml, bulk of 1000 tubes/bag
60-0008	PCR tube, strip of 8, separate cap strip, bulk of 125 strips/bag
60-0012	PCR tube, strip of 12, separate cap strip, bulk of 85 strips/bag


Contact

analysisio GmbH
 Brandteichstrasse 20
 D-17489 Greifswald

Fon: +49(0)3834 550263
 Fax: +49(0)3834 550222
 Mobile: +49(0)171 3406464
 Email: info@analysisio.com



Cardbord freezer boxes



Strong cardboard
 Superior white coating
 Large writing area on the lid
 For upright or chest freezer racks
 Space for 81 or 100 tubes
 Temperature range from -196 °C to 121 °C
 Two available sizes: either 2" or 3" height

Contact


analysisio GmbH
 Brandteichstrasse 20
 D-17489 Greifswald

Fon: +49(0)3834 550263
 Fax: +49(0)3834 550222
 Mobile: +49(0)171 3406464
 Email: info@analysisio.com

Article number	Description
90-1281	2" cryobox with 81 places
90-1200	2" cryobox with 100 places
90-1381	3" cryobox with 81 places
90-1300	3" cryobox with 100 places



Pipet tips



10 μ l, 200 μ l and 1000 μ l pipet tips
 RNase & DNase Free
 Polypropylene
 Natural, yellow and blue color
 Autoclavable
 Graduated
 1000 Tips/Bag
 Low retention and beveled fine point
 Fit to most common research grade pipettor
 Suitable for manual and electronic units

Pipet tip rack




200ul Pipet tips, 96 tips/rack
 RNase & DNase Free
 Polypropylene
 Yellow color
 Autoclavable
 Graduated
 Beveled point

Article number	Description
20-0010	10 μ l natural color pipet tips, bulk pack of 1000/bag
20-0200	200 μ l yellow color pipet tips, bulk pack of 1000/bag
20-1000	1000 μ l blue color pipet tips, bulk pack of 1000/bag
21-0200	200 μ l yellow color pipet tips, rack pack of 96 tips


Contact

analysisio GmbH
 Brandteichstrasse 20
 D-17489 Greifswald

Fon: +49(0)3834 550263
 Fax: +49(0)3834 550222
 Mobile: +49(0)171 3406464
 Email: info@analysisio.com



Transfer pipets



General purpose disposable transfer pipets
 Quick and safe sampling and transfer of fluids
 One piece, flexible
 Polyethylene, 100 % inert
 Suitable for gas sterilization
 Gamma-radiation-sterilized products available
 Graduated or non-graduated
 Many different capacities and shapes (see table)

Article number	Description
30-0604	Non-graduated, V: 5,6 ml, l:148 mm, fine tip: 100 drops/ml , bulk: 500/pack
30-0135	Graduated, V: 6,2 ml, l: 160 mm, 26 drops/ml, bulk: 500/pack
30-0138	Graduated, V: 7,5 ml, l: 160 mm, 22 drops/ml bulk: 500/pack
30-0616	Non-graduated, V: 5,5 ml, l: 232 mm, 28 drops/ml, bulk: 500/pack
Number-A1	Sterile, individually packed in a paper/plastic blister, 500/pack
Number-B20	Sterile, 20 pieces in a bag, 500/pack


Contact

analysisio GmbH
 Brandteichstrasse 20
 D-17489 Greifswald

Fon: +49(0)3834 550263
 Fax: +49(0)3834 550222
 Mobile: +49(0)171 3406464
 Email: info@analysisio.com



Solution basin



Disposable, strong and stable solution basin
 Clear color
 55ml Capacity
 PVC
 Nonsterile
 Disposable
 For single or multi-channel pipetters
 Bulk packed




Sterile solution basin
 Gamma radiation sterilized
 White PS material
 55ml capacity
 Disposable
 For single or multi-channel pipetters
 Individually packed

Contact

analysisio GmbH
 Brandteichstrasse 20
 D-17489 Greifswald

Article number	Description
25-0052	PVC, nonsterile 55 ml solution basin, bulk pack of 100
25-0051	PS, gamma radiation sterilized solution basin, single packed, pack of 50

Fon: +49(0)3834 550263
 Fax: +49(0)3834 550222
 Mobile: +49(0)171 3406464
 Email: info@analysisio.com



Drosophila bottles



High quality PP or PET bottles
 Good Visibility (PET: super clear)
 Molded or printed Graduations
 PP bottles autoclavable
 PET bottles non-autoclavable
 Bulk or tray pack
 Capacity: 177ml
 Dimensions: 57mm l x 57mm w x 103mm h
 Square bottom

Drosophila vials




High quality PS material
 Excellent visibility
 Bulk or tray pack
 Three sizes:
 Wide (28,5 x 95 mm)
 Narrow (25 x 95 mm)
 And short-narrow (25 x 75 mm)
 Non-autoclavable

Article number	Description
51-17720	PP drosophila bottle, bulk of 200/pack, autoclavable
51-17721	PP drosophila bottle, tray with 25 bottles, autoclavable
53-1770	PET drosophila bottle, bulk of 500/pack, super clear, non-autoclavable
53-1771	PET drosophila bottle, tray with 25 bottles, super clear, non-autoclavable
51-0500	PS drosophila vial, narrow, bulk of 500/pack
51-0501	PS drosophila vial, narrow, tray of 100 vials
51-0700	PS drosophila vial, short-narrow, bulk of 500/pack
51-0701	PS drosophila vial, short-narrow, tray of 100 vials
51-0800	PS drosophila vial, wide, bulk of 500/pack
51-0801	PS drosophila vial, wide, tray of 100 vials

Contact

analysisio GmbH
 Brandteichstrasse 20
 D-17489 Greifswald

Fon: +49(0)3834 550263
 Fax: +49(0)3834 550222
 Mobile: +49(0)171 3406464
 Email: info@analysisio.com



Tri-cornered beakers



Unbreakable, chemical inert multipurpose beakers
 6 different volumes with molded graduations
 Three pour spots
 Disposable as well as reusable
 Autoclavable

Petri dishes



High quality PS material
 Designed for excellent visibility
 Sterile
 Triple-vented lid for less condensation
 Stackable for your convenience
 Dimensions: 100mm x 15mm

Article number	Description
93-0050	Tri-cornered beaker, 50 ml, pack of 25 in a sleeve
93-0100	Tri-cornered beaker, 100 ml, pack of 25 in a sleeve
93-0250	Tri-cornered beaker, 250 ml, pack of 25 in a sleeve
93-0400	Tri-cornered beaker, 400 ml, pack of 25 in a sleeve
93-0800	Tri-cornered beaker, 800 ml, pack of 25 in a sleeve
93-1000	Tri-cornered beaker, 1000 ml, pack of 25 in a sleeve
66-1501	Petri-dishes, 100 mm x 15 mm, sterile, pack of 20 in a sleeve

Contact

analysisio GmbH
 Brandteichstrasse 20
 D-17489 Greifswald

Fon: +49(0)3834 550263
 Fax: +49(0)3834 550222
 Mobile: +49(0)171 3406464
 Email: info@analysisio.com



Cell scrapers and lifters



Disposable PE cell scrapers and lifters
 Molded handle in two different lengths
 Pivoting blade
 Gamma-radiation sterilized
 Non-pyrogenic
 Individually wrapped

Embedding cassettes



Embedding cassettes and spare O-rings
 For biopsy or tissue processing
 High density polymer
 Efficient flow-through strip or square holes
 Resistant to chemical action of histological solvents
 45° angle writing area
 Different colors

Article number	Description
70-1180	Cell scraper, l: 180 mm, 18 mm blade, sterile, bulk of 100 pieces
70-1250	Cell scraper, l: 250 mm, 18 mm blade, sterile, bulk of 100 pieces
70-2180	Cell lifter, l: 180 mm, 20 mm blade, sterile, bulk of 100 pieces
41-100X	Embedding cassettes, strip holes for tissue processing, bulk of 500 pieces
41-200X	Embedding cassettes, square holes for biopsy processing, bulk of 500 pieces
41-300X	Embedding spare O-rings, bulk of 500 pieces
Color code: X=1: white, X=2: pink, X=3: green, X=4: yellow, X=5: blue	

Contact

analysisio GmbH
 Brandteichstrasse 20
 D-17489 Greifswald

Fon: +49(0)3834 550263
 Fax: +49(0)3834 550222
 Mobile: +49(0)171 3406464
 Email: info@analysisio.com

Slide storage boxes



Strong, durable plastic
 Cork or foam lining
 Superior design protects against dust and damage
 Space for 25 or 100 3 X 1 inch slides
 Assorted colors
 Stackable

Slide mailers



Strong, durable cardboard slide mailer / holder
 Safe slide storage and transport
 Space for 1, 2 or 20 slides
 20 slides version:
 Premium glossy finish
 Printed index
 Sturdy, efficient storage box for 10 holders

Article number	Description
41-5025	25 slides storage box, foam lining, 5,5 X 3,75 X 1,25 inch, pack of 5 pieces
41-5100	100 slides storage box, foam lining, 8,25 X 6,37 X 1,25 inch, pack of 5 pieces
41-6025	25 slides storage box, cork lining, 5,5 X 3,75 X 1,25 inch, pack of 5 pieces
41-6100	100 slides storage box, cork lining, 8,25 X 6,37 X 1,25 inch, pack of 5 pieces
41-4001	Mailer for 1 slide, pack of 100
41-4002	Mailer for 2 slides, pack of 100
41-4020	Mailer / holder for 20 slides, box of 10

Contact

analysisio GmbH
 Brandteichstrasse 20
 D-17489 Greifswald

Fon: +49(0)3834 550263
 Fax: +49(0)3834 550222
 Mobile: +49(0)171 3406464
 Email: info@analysisio.com

Contact and Ordering Information

Please feel so free to contact us if you have any question or if you need more Information.

We like to help you with more information about BIOLOGIX laboratory supplies.

To contact us you may use the contact form on our homepage www.analysisio.com or write us an email to info@analysisio.com.

To request a quotation or to place an order you may use the order/quotation form on our homepage www.analysisio.com or write us an email to info@analysisio.com.

Please use the article numbers of this catalogue for ordering or requesting a quotation.

We look forward to your order and would like to send you an offer.

Of course you may contact us also via phone or fax.

You find an actual pricelist in the download area where you found this catalogue.

Imprint

analysisio GmbH

Brandteichstrasse 20

D-17489 Greifswald, Germany

Tel.: +49 (0)3834 550263

Fax: +49 (0)3834 550222

E-mail: info@analysisio.com

General manager: Jens Pörksen

Commercial registry: Stralsund HRB 7139

Sales tax id. number: DE252491041

