

analysio® Reference Electrode Container

More security and comfort
Less handling and maintenance
No contamination of electrode and sample



analysio® Reference Electrode Container

The analysio® Reference Electrode Container avoids as well contamination of the reference electrode as of the sample.

A gel electrolyte bridge connects the reference electrode reliably with the sample.


The choice of the suitable electrolyte and the easy change of the gel electrolyte bridge ensure interference and trouble free measurements.

The reference electrode is always stored under optimal conditions. The electrolyte in the analysio® Reference Electrode Container is easily fitted to the used reference electrode.

The reference electrode does not need more space and handling than the gel electrolyte bridge.

analysio® Reference Electrode Container		
Type	Ø reference electrode	Filling electrolyte
REC 012-01	12 mm	KCl
REC 004-02	4 mm	K ₂ SO ₄
REC 012-01	12 mm	KCl
REC 004-02	4 mm	K ₂ SO ₄

Gel electrolyte bridges		
Type	Electrolyte (c = 1,0 mol/l)	Anwendungsbeispiel
REC 000-01	KCl	Standard
REC 000-02	NaCl	Physiology*
REC 000-03	KNO ₃	Measurements with ion selektive electrodes
REC 000-04	NaNO ₃	
REC 000-05	KHPO ₄	
REC 000-06	NaHPO ₄	



Measuring with the analysio® Reference Elektrode Container

Contact

analysio GmbH
 Brandteichstrasse 20
 D-17489 Greifswald

Fon: +49(0)3834 550263
 Fax: +49(0)3834 550222
 Mobile: +49(0)171 3406464
 Email: info@analysio.com



Contact and Ordering Information

Please feel so free to contact us if you have any question or if you need more Information.

We like to help you with more information about the analysisio® Reference Electrode Container.

To contact us you may use the contact form on our homepage www.analysisio.com or write us an email to info@analysisio.com.

To request a quotation or to place an order you may use the order/quotation form on our homepage www.analysisio.com or write us an email to info@analysisio.com.

Please use the article numbers of this catalogue for ordering or requesting a quotation.

We look forward to your order and would like to send you an offer.

Of course you may contact us also via phone or fax.

You find an actual pricelist in the download area where you found this catalogue.

Imprint

analysisio GmbH
 Brandteichstrasse 20
 D-17489 Greifswald, Germany
 Tel.: +49 (0)3834 550263
 Fax: +49 (0)3834 550222
 E-mail: info@analysisio.com
 General manager: Jens Pörksen
 Commercial registry: Stralsund HRB 7139
 Sales tax id. number: DE252491041

Contact

analysisio GmbH
 Brandteichstrasse 20
 D-17489 Greifswald

Fon: +49(0)3834 550263
 Fax: +49(0)3834 550222
 Mobile: +49(0)171 3406464
 Email: info@analysisio.com

